


Article

Monitoring the Chloride Concentration in International Scheldt River Basin District Water Using a Low-Cost Multifunction Data Acquisition Board

Wanda J. Guedens ^{1,2,*} , Monique Reynders ² , Koen Van Vinckenroye ², Jan Yperman ^{2,3}  and Robert Carleer ^{2,3}

¹ Biomolecule Design Group, Institute for Materials Research, Hasselt University, Campus Diepenbeek, Agoralaan-Building D, BE-3590 Diepenbeek, Belgium

² Faculty of Sciences, Hasselt University, Campus Diepenbeek, Agoralaan-Building D, BE-3590 Diepenbeek, Belgium; monique.reynders@uhasselt.be (M.R.); koen.vanvinckenroye@uhasselt.be (K.V.V.); jan.yperman@uhasselt.be (J.Y.); robert.carleer@uhasselt.be (R.C.)

³ Applied and Analytical Chemistry, Institute for Materials Research, Hasselt University, Campus Diepenbeek, Agoralaan-Building D, BE-3590 Diepenbeek, Belgium

* Correspondence: wanda.guedens@uhasselt.be; Tel.: +32-11-268324

Received: 20 June 2018; Accepted: 30 July 2018; Published: 2 August 2018



Abstract: In analytical chemistry laboratories, to gather in the shortest time as many data as possible with the utmost accuracy and precision, high throughput automated setups are indispensable. In the present study, to determine the chloride concentration in the international Scheldt river basin district, experiments are carried out utilizing a thermostatically controlled semi-automated setup. A novel ICT-based method is developed using a low-cost multifunction Data Acquisition Board (DAQ) controlled by a homebuilt LabVIEW™ program. Specifically, this approach enables a correlation between different parameters i.e., droplet volume, temperature, A/D voltage conversions. Here, processing experimental data of a potentiometric precipitation titration utilizing a silver nitrate standard solution as titrant in a manual burette equipped with a controllable electronic valve allows for a preliminary indication of the titration end point via the Virtual Instrument (VI) numerical first derivative tool in the LabVIEW software. The LabVIEW tool is compared with the well-known Gran method implemented in the LabVIEW program, emphasizing an accurate performance of the setup to determine the chloride concentration in fresh river water. We are confident that our findings are evidence of the versatile and powerful features of the LabVIEW controlled DAQ in the analytical chemistry laboratory.

Keywords: electroanalytical chemistry; potentiometric titration; ionic solutions; multifunction data acquisition; chloride; solubility product; natural surface water; fresh river water

1. Introduction

According to the Council Directive of the European Communities concerning the optimal quality for natural surface water, the average threshold value for the chloride concentration is 5.642 mmol/L (200 mg/L) [1,2]. In natural surface waters, the chloride concentration depends not only on geological and geographical background conditions, but on ecological circumstances as well [3]. On the one hand, high chloride concentrations in water are not known to have toxic effects on humans, although large amounts may act corrosively and be harmful to plant life. Chlorides are present in both fresh and salt

water and, are essential elements of life [4]. On the other hand, salinization of fresh water and dry periods pose a serious problem for the surrounding areas of river and stream basins.

A high chloride concentration may be due to environmental contamination. In the present study, the experimental settings are validated to consolidate the determined chloride concentrations in natural surface water [5].

In our analytical chemistry labs, the accuracy of low-cost automated potentiometric measurements e.g., determination of common salt concentration, has increased over time [6]. In fundamental and applied chemistry research in general, it is of the utmost importance to record and process physicochemical parameters e.g., the pH of a buffer solution, with the highest accuracy and precision. For these measurements, the concentration of hydrogen ions was measured by a hydrogen ion selective electrode [7]. These results encouraged us to measure chloride ions in fresh water using a Silver/Sulfide Solid State ion selective electrode in an adapted potentiometric titration setup.

A powerful semi-automated analytical setup was developed for monitoring multiple variables using a low-cost multifunction Data Acquisition Board (DAQ) [8] and a homemade LabVIEW™ program. In this study, the chloride concentration in a temperature-controlled, unknown NaCl solution, as well as in natural surface water, is determined. On a hardware level, the semi-automated setup is equipped with a temperature sensor and an electronic valve controlled by the DAQ, while a pH/mV potentiometer and a cryostat to maintain a pre-set temperature are connected with the computer via a RS 232 interface port. On the software level, in this laboratory environment, the potential, $E(\text{Ag}^+)$, is monitored during a potentiometric precipitation titration of a primary standard NaCl solution with a standard AgNO_3 solution using a pH/mV potentiometer. At the end point of the titration, $E(\text{Ag}^+)$ and the added volume of the AgNO_3 ($V(\text{Ag}^+)$) are determined by the Virtual Instrument (VI) numerical first derivative tool in the LabVIEW™ 2015 SP 1 (version 15.0.1) software. Via the well-known titration formula ($c(\text{Cl}^-)$) can be calculated. Via the calibration curve, $E(\text{Ag}^+)$ versus $-\log(c(\text{Ag}^+))$ for different thermostatically controlled AgNO_3 standard solutions, the equilibrium concentration of Ag^+ i.e., $[\text{Ag}^+]$ at the end point is derived, and thus, the solubility product of silver chloride (K_{sp}) at 298 K.

In view of an increase of accuracy in determining the endpoint in the precipitation titration curve by the Gran plot mathematical approach, a homemade LabVIEW program is developed.

Thus, the envisaged Cl^- concentration ($c(\text{Cl}^-)$) in standard conditions and in fresh river water is indirectly calculated from the measured $E(\text{Ag}^+)$ data as a function of $V(\text{Ag}^+)$. Note that for the Gran plot equations [9], the above-mentioned calibration curve is still crucial to convert the measured $E(\text{Ag}^+)$ values into $[\text{Ag}^+]$. The Gran method is validated by comparing the independently obtained analytical concentrations of silver and chloride ions and indirectly K_{sp} , via the LabVIEW VI tool, based on a single measuring point at the end, and the Gran plot strategy for a range of experimental titration data near the end point. This not only confirms the reliability of the tailor-made Gran plot program, but also the reproducibility of the experimental conditions of the semi-automated setup for the potentiometric determination of the chloride concentration in fresh river water [10].

To summarize, our method enables the simultaneous determination of the free/analytical concentrations of chloride and silver ions, and the solubility product at a given temperature, in one single experiment.

It is shown that the values of the experimentally determined chloride concentration and the related solubility product of silver chloride by the Gran method are more reproducible and consistent for a series of titrations than by the VI method, thus emphasizing an accurate performance of the used semi-automated potentiometric setup for the analysis of fresh river water.

2. Materials and Methods

2.1. Chemicals

Forty milliliters natural surface water of the Demer river (50°55'28.7" N, 5°23'22.8" E), a tributary of the international Scheldt river basin district, was sampled using a dedicated vessel [11,12].

Solid silver nitrate, AnalaR Normapur[®], sodium chloride, AnalaR Normapur[®] and, Ionic Strength Adjuster solution (ISA, 5 M NaNO₃) for silver measurements, Cat. No. 940011, Thermo Scientific were supplied by VWR International, Oud-Heverlee, Belgium. Water was purified by an arium[®] 611DI | 611UV apparatus (VWR International, Oud-Heverlee, Belgium).

2.2. Equipment and General Setup

In the framework of the presented study, the traditional manual configuration of a titration is customized to a semi-automated setup (Figure 1 and Figure S1).




Figure 1. Configuration of the semi-automated setup for thermostatically controlled potentiometric measurements.

A homemade double-walled titration vessel insulated with a polyurethane jacket (Figure 2) was connected to a cryostat (Figure S2) (Polyscience Digital Temperature Controller (model 9110), VWR International, Oud-Heverlee, Belgium).



Figure 2. A homemade double walled titration vessel (left); a polyurethane jacket (right).

The computer was extended with a multifunction NI PCI-6220-DAQ board, a homemade power supply interface box and the National Instruments Software Platform Bundle Spring 2016 (NI corp.) [8]. A temperature sensor, PT1000 (Schott type W 5780 NN supplied by VWR International, Oud-Heverlee,

Belgium) to measure the temperature of the solution in the vessel (Figure S3) was connected via DAQ built into a HP Compaq Elite 8300 Micro tower working with the operating system Microsoft Windows 10. The temperature of the water in the cryostat, connected via a RS 232 interface port, and thus, also of solution in the insulated titration vessel, was automatically set to 298 K. The temperature was monitored to allow the cryostat to respond to the least temperature change.

The sodium chloride solution in the first part of the experiment, and the filtered surface water in the second part, were added with a 20 mL (± 0.03 mL, 20 °C) bulb pipette, BLAUBRAND® to the insulated titration vessel placed on a magnetic stirrer (type IKA Mini-MR standard 100–240 V, 50/60 Hz, 0–1500 rpm, supplied by VWR International, Oud-Heverlee, Belgium). The surface water was filtered on GN-6 Metrical® S-Pack membrane Disc Filters Prod. No. 66278, using Sentino™ Magnetic Filter Funnels, Part Number 4271, PALL® Life Sciences to remove insoluble organic and inorganic particles. The solution in the insulated titration vessel was constantly stirred with a magnetic stirring rod at a low velocity to minimize the capture of carbon dioxide from the air.

A 50 mL (± 0.05 mL, 20 °C) burette BLAUBRAND® filled with the silver nitrate solution was manually assembled just above the stopcock using inert silicone glue-elastosil (supplied by VWR International, Oud-Heverlee, Belgium) to an electronic valve, 2-way solenoid valve, normally open, STV-02-M6(1/4U)KG-12 V DC, Metrohm, Antwerpen, Belgium. A Silver/Sulfide Solid State Half-Cell Ion Selective Electrode (ISE) Orion 9416SC Thermo Scientific [13] (supplied by Thermo Fisher Scientific, Geel, Belgium) and a Ceramic Hg/Hg₂SO₄ reference electrode (+15 °C . . . +40 °C), SI Analytics B 3610+ (supplied by VWR International, Oud-Heverlee, Belgium), introduced in the vessel, are connected to a potentiometer pH/mV meter WTW inoLab® pH 730 BNC-1A20-210 (supplied by VWR International, Oud-Heverlee, Belgium). The voltage measurements of the potentiometer were converted by the DAQ, controlled by a homemade LabVIEW™ program.

2.3. Volume Control of Titrant Reagent in the Semi-Automated Setup

Volume control of the titrant reagent is crucial in the semi-automated setup, since the measurements during the titration are executed in function of the added volume. Therefore, a 50 mL burette was equipped with a controlled electronic valve, as described in the general setup. The isolated solenoid diaphragm valve was chosen for its inertness to most chemicals, alkaline as well as acid. Highly inert material, polytetrafluoroethylene (PTFE), was used for the parts in contact with the fluid, resulting in a design that can be used for a wide range of fluids.

The solenoid diaphragm valve has a function 'Normal Open' to fill the burette with the titrant reagent. To close the valve, the controlling electronic power supply circuit of the homemade interface box (Figure 3) provides 12 V DC across the solenoid (COIL).




Figure 3. Electronic power supply circuit of the homemade interface box.

To calibrate the volume of the droplets dripping from the burette, a LabVIEW program was written to automatically drip a predetermined number of droplets of the solution into the vessel in a set time frame (3 droplets/0.8 s). The volume was manually read from the burette and registered into the program. This procedure needs to be repeated at the beginning of each experiment.

2.4. Automated Registration of the Ag^+ -ion Activity

For the automated registration of the Ag^+ ion activity at the Ag/S^{2-} electrode with respect to the $\text{Hg}/\text{Hg}_2\text{SO}_4$ electrode, voltage (mV) data were read out by the computer via RS232 communication with the potentiometric mV meter.

2.5. The Homemade LabVIEW Program

Currently, in our lab, the LabVIEW run time engine 2013 [8] is installed as a communication tool between the homemade LabVIEW program and the DAQ. The homemade program using LabVIEW[®] software (supplied by National Instruments Belgium nv, Zaventem, Belgium), combines different tasks: Controlling the valve, calculating the added volume via calibration volume droplets of the titrant, reading the potentiometric mV meter and data processing.

2.6. Methodology

In Scheme 1, the methodology is outlined. In **Part I** of the experimental procedure, a calibration of the droplet volume was performed by the homemade LabVIEW program controlling the normally open valve connected to the burette. The valve opens and closes automatically five times while the AgNO_3 solution is dripping out. The volume poured out was read from the burette, and the program calculated the mean droplet volume from a triplicate measurement. In **Part II** of the experimental procedure, the linear regression line (Nernst equation) $E(\text{Ag}^+)$ versus $-\log c(\text{Ag}^+)$ was constructed via the modelling LabVIEW tool, starting from the mV-measurements of five different, thermostatically controlled (298 K) standard AgNO_3 solutions adjusted with ISA (1 volume ISA/50 volumes ionic solution) (Figure 4). The Ag^+ activity of a series of AgNO_3 solutions in the homemade double walled vessel was measured via the Ag/S^{2-} electrode.





Figure 4. Linear regression line, E versus $-\log c(\text{Ag}^+)$ for the standardization of a silver nitrate solution.

Subsequently, in **Part III** of the experimental procedure, 40 mL of a thermostatically controlled 0.005 M NaCl solution adjusted with ISA was titrated with an unknown 0.01 M AgNO_3 solution (Scheme 1 **Part III a**). Finally, 40 mL of a filtered, thermostatically controlled sample of natural surface

river water adjusted with ISA solution was titrated with the standardized 0.01 M AgNO_3 solution while stirring at a constant rate (Scheme 1 **Part III b**).

The registered titration curve (E versus $V(\text{AgNO}_3)$) was processed by the home made LabView program via the VI method (**Part IV**) and the Gran plot method (**Part V**). Particularly, in **Part IV**, the VI tool was used to determine the different variables, $V(\text{Ag}^+)$ and $[\text{Ag}^+]$ at endpoint (EP), $c(\text{Cl}^-)$ (via titration formula), and K_{sp} . In **Part V**, the former variables were calculated via the Gran equations and accompanying Gran plot.



Scheme 1. Flow chart of the experimental procedure.

3. Results & Discussion

The experimental setup for a manual precipitation titration was modified to create a semi-automated setup by using a computer configuration enhanced with a DAQ. A homemade LabVIEW program controlled the hyphenated peripherals i.e., the electronic valve for volume measurements, the temperature-dependent resistance sensor, a potentiometric pH/mV meter and the cryostat for maintaining a constant temperature in the homemade vessel. Hereby, different experimental conditions and variables are optimized simultaneously. In this improved environment, the analytical chloride concentration was determined in a sample of fresh river water via a potentiometric silver chloride precipitation titration as described above.

Subsequently, a standard silver nitrate solution was used as a titrant for the determination of the analytical chloride concentration in an unknown NaCl solution, and finally in fresh river water.

The titration curves, E (mV) versus $V(\text{AgNO}_3)$ (mL), of the titration of an unknown NaCl solution and fresh river water, respectively (Figures 5a and 6a) provide information on the end point derived by the numerical first derivative method, the VI tool within the LabVIEW software, as well as by the Gran plot equations in the tailor-made LabVIEW program (Tables 1 and 2). The VI method is limited due to the numerical character of the derivatives leading to one single titration end point: The slightest deviation in the potential, E , at the end point of the titration has a dramatic impact on the value of K_{sp} . The Gran plot strategy is not subject to this drawback because of the use of a range of experimental titration data near by the end point.

Table 1. Results at the end point of a precipitation titration of three replicates of a 40.00 ± 0.04 mL sample at 298 K of a 0.005 M NaCl solution with a 0.01 M standard AgNO_3 solution by using the LabVIEW Virtual Instrument (VI) method versus the Gran method for the added volume (mL) AgNO_3 titrant; E (mV) derived from the maximal numerical derivative of the titration curve; the solubility product of silver chloride (K_{sp}) derived from the remaining free concentration of Ag^+ ions and Cl^- ions using the LabVIEW VI method versus the Gran method [14].


Experiment	T (K) at the End Point	V(Ag^+) (mL)	LabVIEW VI Method		Gran Method			
			E ** (mV)	$K_{\text{sp}} \times 10^{10}$	E ** (mV)	c(Ag^+) * (mmol/L)	c(Cl^-) * (mmol/L)	$K_{\text{sp}} \times 10^{10}$
Titration 1	298.45	20.16	−123.8	1.19	−114.9	8.97	4.52	2.41
Titration 2	297.81	20.02	−129.5	0.76	−115.1	8.62	4.31	2.37
Titration 3	298.24	20.31	−124.2	1.15	−115.3	8.54	4.33	2.33
Mean				1.03		8.71	4.39	2.37
Standard deviation (St.dev.)				0.24		0.23	0.11	0.04

Note: * Analytical concentration. ** E value at the end point of the titration.


Table 2. Results at the end point of a precipitation titration of three replicates of a 40.00 ± 0.04 mL sample of fresh river water at 298 K by using the LabVIEW VI method versus the Gran method.

Experiment	T (K) at the End Point	V(Ag^+) (mL)	LabVIEW VI Method		Gran Method			
			E ** (mV)	$K_{\text{sp}} \times 10^{10}$	E ** (mV)	c(Ag^+) * (mmol/L)	c(Cl^-) * (mmol/L)	$K_{\text{sp}} \times 10^{10}$
Sample 1	298.01	6.42	−134.4	0.517	−118.0	6.65	1.07	1.88
Sample 2	298.29	6.20	−122.5	1.32	−117.8	6.42	1.00	1.91
Sample 3	298.31	6.57	−124.4	1.14	−117.1	6.34	1.04	2.02
Mean				0.99		6.47	1.04	1.94
St.dev.				0.42		0.16	0.04	0.07

Note: * Analytical concentration [2]. ** E value at the end point.





(a)



(b)

Figure 5. (a) Precipitation titration curve (titration 2; Table 1) of an unknown NaCl solution at 297.81 K (temperature at the end point), E (mV) versus V (mL) AgNO₃ standard solution (red); plot of numerical first derivatives (built-in LabVIEW VI function) with a maximum at the endpoint (orange); (b) Gran plot with an intersection of two linear regression lines at the end point of the titration.



(b)

Figure 6. (a) Precipitation titration curve (titration 2; Table 2) of a sample of fresh river water at 298.29 K, E (mV) versus V (mL) standard AgNO₃ solution (red); plot of numerical first derivatives (built-in LabVIEW VI function) with a maximum at the endpoint (orange); (b) Gran plot with an intersection of two linear regression lines at the end point of the titration.

Indeed, it is clear from Tables 1 and 2 that the variation in the individual K_{sp} values derived from the Gran plot is less pronounced compared to the one from the VI method. Additionally, from the Gran plot beyond the EP, the analytical concentration of silver ions, as well as that of chloride ions, is derived.

In the Supplementary Materials, an Excel (2016 MS Office) file is provided with different worksheets i.e., 'Calibration curve', 'Standardization of NaCl', 'Titration of fresh river water samples'. The Gran plot equations are used to calculate the different variables i.e., the analytical concentration of Ag^+ and Cl^- and K_{sp} in Tables 1 and 2.

In the Gran plot, the two intersecting linear regression lines before and beyond the EP are derived as follows [15]



At equilibrium

$$K_{sp} = [Ag^+] \cdot [Cl^-] \quad (2)$$

The equilibrium concentration of the silver ions at the EP, $[Ag^+]$ was derived from the calibration curve (Figure 4):

$$E = a \cdot (-\log [Ag^+]) + b$$

$a = \text{slope}$; $b = \text{the intercept}$

$$\frac{1}{[Ag^+]} = 10^{(E-b)/a} \text{ or } [Ag^+] = 10^{(b-E)/a} \quad (3)$$

At the EP the Ag^+ and the Cl^- concentration are equal. K_{sp} was determined by the square of the concentration $[Ag^+]$ and, thus:

$$K_{sp} = [Ag^+]^2 \text{ or } 10^{2(b-E)/a} \quad (4)$$

Derivation of the linear regression line prior to the EP:

$$[Cl^-] = \frac{c(Cl^-) \cdot V - c(Ag^+) \cdot v}{V + v} \quad (5)$$

with $V = \text{initial volume of } Cl^-$ and $v = \text{added volume of } Ag^+$.

Insert Equation (5) into Equation (2):

$$\frac{K_{sp}}{[Ag^+]} = \frac{c(Cl^-) \cdot V - c(Ag^+) \cdot v}{V + v} \quad (6)$$

Considering Equation (3) into Equation (6):

$$(V + v) \cdot 10^{(E-b)/a} = \frac{c(Cl^-) \cdot V - c(Ag^+) \cdot v}{K_{sp}} \quad (7)$$

thus:

$$(V + v) \cdot 10^{E/a} = \frac{-c(Ag^+) \cdot 10^{b/a}}{K_{sp}} \cdot v + \frac{c(Cl^-) \cdot 10^{b/a}}{K_{sp}} \cdot V \quad (8)$$

In Equation (8), the slope of the linear curve,

$$(V + v) \cdot 10^{E/a} \text{ versus } v, \text{ equals } \frac{-c(Ag^+) \cdot 10^{b/a}}{K_{sp}} \quad (9)$$

Derivation of the linear regression line beyond the EP:

$$[Ag^+] = \frac{c(Ag^+) \cdot v - c(Cl^-) \cdot V}{V + v} \quad (10)$$

Considering Equation (3):

$$(V + v) \cdot 10^{-(E-b)/a} = c(\text{Ag}^+) \cdot v - c(\text{Cl}^-) \cdot V \quad (11)$$

thus:

$$(V + v) \cdot 10^{-E/a} = c(\text{Ag}^+) \cdot 10^{-b/a} \cdot v - c(\text{Cl}^-) \cdot 10^{-b/a} \cdot V \quad (12)$$

In Equation (12), the slope of the linear curve

$$(V + v)10^{-E/a} = c(\text{Ag}^+) \cdot 10^{-b/a} \quad (13)$$

It should be noted that no grid data are ignored in the neighborhood of the equivalence point, while at the beginning and the end of the titration curve, 60% of the data are omitted to obtain two intersecting regression lines (Figures 5b and 6b) with optimal correlation coefficients [16]. K_{sp} is calculated by Equation (8). Additionally, the analytical concentration of silver ions, as well as the chloride concentration, is calculated via the slope and intercept of Equation (12), respectively. Alternatively, the chloride concentration can be calculated by means of the volume titrant at the intersection of the two linear regression lines in the Gran plot.

The values of K_{sp} calculated via the LabVIEW VI method and those calculated via the linear regression line prior to the EP in the Gran plot show significant differences (Table 1). However, the values in far-right column (Table 1) are more consistent (lower standard deviation) in the three experiments. Based on these observations, the Gran method is preferred for the quantitative determination of the chloride concentration in fresh river water (Table 2). This is evidenced by the consistency of the K_{sp} values in the rightmost column. Additionally, due to these findings, the analytical concentration of the titrant ($c(\text{Ag}^+)$) has been determined with great certainty via the linear regression line beyond the EP in the Gran plot. Three constraints in the laboratory test conditions were taken into account: The accuracy of the potentiometer (± 0.3 mV), the limitation of the Ag/S^{2-} working electrode to measure at very low $[\text{Ag}^+]_{EP}$ ($\pm 1 \times 10^{-5}$ M), and the ambient temperature surrounding the semi-automated setup that could influence the measuring results. Moreover, the assessment of the reliability of the LabVIEW VI method versus the Gran method revealed that since the exact (bio)chemical composition of fresh river water was unknown, the presence of e.g., possible interfering ions could be ignored.

Consequently, the calculated analytical chloride concentration in fresh river water is acceptable with good certainty.

4. Conclusions

In the present study, a semi-automated setup for potentiometric precipitation titrations was proven to be a valuable analytical tool in ambient laboratory conditions. The multiple characteristics of a series of data can be monitored in real time thanks to the versatile features of the low-cost multifunction DAQ in an enhanced computer configuration. The presented method is a unique combination of a traditional potentiometric setup design for a precipitation titration of chloride ions with a silver ion solution and an advanced automated setup. Moreover, our novel method enables the simultaneous determination of a variety of variables in one single experiment.

Here, the chloride concentration in an unknown NaCl solution and fresh river water is determined to benchmark the numerical VI tool integrated in the LabVIEW software with a homemade LABVIEW Gran plot program. The values of K_{sp} calculated via the LabVIEW VI method and those calculated via the Gran method show significant differences. Based on these observations, the Gran method is preferable for the quantitative determination of the chloride concentration in fresh river water. Additionally, the analytical concentration of the titrant ($c(\text{Ag}^+)$) was determined with good certainty.

The chloride concentration was found to be far below the threshold value, and thus, in accordance to the Council Directive of the European Communities concerning the optimal quality for natural

surface water. In other words, the Scheldt river basin district in the area of the Demer river still does not suffer from salinization.

Conclusively, the procedure benefits from the versatility of the DAQ as a controlling and monitoring low-cost device, and the assets of the homemade LabVIEW program.

The presented semi-automated setup paves the way for the use of controlling peripherals to acquire series of data, as well as for homemade programming in LabVIEW.

Supplementary Materials: The following are available online at <http://www.mdpi.com/2073-4441/10/8/1025/s1>, Figure S1: Computer configuration, Figure S2: A cryostat, type Polyscience Digital Temperature Controller (model 9110), Figure S3: Electronic diagram of the PT1000, Schott type W 5780 NN (**top**); the A/D converter on the NI PCI-6220 DAQ Board (**bottom**), MS Office Excel S1: OVERVIEW_TITRATIONS_FINAL_Guedens & Reynders.

Author Contributions: Conceptualization, W.J.G. and M.R.; data curation, W.J.G. and M.R.; methodology, W.J.G. and M.R.; software, K.V.V. and J.Y.; validation, J.Y. and R.C.; formal analysis, W.J.G. and M.R.; investigation, W.J.G. and M.R.; writing-original draft preparation, W.J.G. and M.R.; writing-review & editing, R.C.; visualization, W.J.G., M.R. and K.V.V.

Funding: This research received no external funding.

Conflicts of Interest: The authors declare no conflict of interest.

References

1. ANNEX 3 to the Commission Staff Working Document Accompanying the Report from the Commission in Accordance with Article 3.7 of the Groundwater Directive 2006/118/EC on the Establishment of Groundwater Threshold Values—Information on the Groundwater Threshold Values of the Member States. Available online: http://ec.europa.eu/environment/water/water-framework/groundwater/pdf/com_swd_annex_iii.pdf (accessed on 16 May 2018).
2. Chloride in Drinking-Water Background Document for Development WHO Guidelines for Drinking-Water Quality. Available online: www.who.int/water_sanitation_health/dwq/chloride.pdf (accessed on 23 May 2018).
3. Walton, R.K.; Samuel, V.P.; Hackley, K. *The Sources, Distribution, and Trends of Chloride in the Waters of Illinois*; Bulletin B-74; Illinois State Water Survey: Champaign, IL, USA, 2012.
4. Shukla, M.; Arya, S. Determination of chloride ion (Cl^-) concentration in Ganga river water by Mohr Method at Kanpur, India. *Green Chem. Technol. Lett.* **2018**, *4*, 6–8. [[CrossRef](#)]
5. Jakmune, J.; Patimapornlert, L.; Suteerapataranon, S.; Lenghor, N.; Grudpan, K. Sequential injection with lab-at-valve (LAV) approach for potentiometric determination of chloride. *Talanta* **2004**, *65*, 789–793. [[CrossRef](#)] [[PubMed](#)]
6. Yperman, J.; Carleer, R.; Reggers, G.; Mullens, J.; Van Poucke, L. Automation of potentiometric measurements—Determination of water-extractable sodium in bread using a sodium-ion selective electrode with minimum sample preparation. *J. AOAC Int.* **1993**, *76*, 1138–1148.
7. Velghe, I.; Carleer, R.; Yperman, J.; Schreurs, S.; D’Haen, J. Characterization of adsorbents prepared by pyrolysis of sludge and sludge/disposal filter cake mix. *Water Res.* **2012**, *46*, 2783–2794. [[CrossRef](#)] [[PubMed](#)]
8. National Instruments. Available online: www.ni.com (accessed on 14 March 2018).
9. Gans, P.; O’Sullivan, B. GLEE, a new computer program for glass electrode Calibration. *Talanta* **2000**, *51*, 33–37. [[CrossRef](#)]
10. Jeffery, G.H.; Bassett, J.; Mendham, J.; Denney, R.C. *Vogel’s Textbook of Quantitative Chemical Analysis*, 5th ed.; Longman Scientific & Technical Longman Group: Harlow Essex, UK, 1989; pp. 340–359, ISBN 0-470-2151.
11. International River Basin District of the Scheldt. Available online: http://carto1.wallonie.be/webgis_escaut_public_en/viewer.htm (accessed on 13 June 2018).
12. Baeyens, W.; van Eck, B.; Lambert, C.; Wollast, R.; Goeyens, L. General description of the Scheldt estuary. *Hydrobiologia* **1998**, *366*, 1–14. [[CrossRef](#)]
13. Thermo Fisher Scientific Inc. *Thermo Scientific User Guide Silver/Sulfide Ion Selective Electrode*; Thermo Fisher Scientific Inc.: Waltham, MA, USA, 2009.
14. Berger, M. Potentiometric Determination of Chloride in Natural Waters: An Extended Analysis. *J. Chem. Educ.* **2012**, *89*, 812–813. [[CrossRef](#)]

15. Location of End Point in Potentiometric Argentometric Titration Using Gran Plot and Titration Errors. Available online: <https://www.scribd.com/doc/38516302/LOCATION-OF-END-POINT-IN-POTENTIOMETRIC-ARGENTOMETRIC-TITRATION-USING-GRAN-PLOT-AND-TITRATION-ERRORS> (accessed on 14 March 2018).
16. Potentiometry in Analytical Chemistry. Available online: <http://www.nepachemistry.com/2011/03/potentiometry-in-analytical-chemistry.html> (accessed on 14 March 2018).



© 2018 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (<http://creativecommons.org/licenses/by/4.0/>).