


Towards lab-on-card, label-free immunoassays based on electrical impedance spectroscopy

A high throughput screening device

Belgian Physical Society General Scientific Meeting,
ULBruxelles, 22 May 2019

Frederik Vreys



Ministerium für Wirtschaft, Energie,
Industrie, Mittelstand und Handwerk
des Landes Nordrhein-Westfalen



Provincie Noord-Brabant



Ministerie van Economische Zaken



AGENTSCHAP
INNOVEREN &
ONDERNEMEN

Vlaanderen
is ondernemen

HEALTHY AGING




Contents







- Introduction
 - Healthy aging
 - Sensor principle
- Characterization
 - Contact angle measurements
 - Antigen coating on substrates
- Sensor development
 - Lab-on-card principle
 - Lab-on-card with impedimetric cell
 - Lab-on-card with microplate format
 - Multichannel impedimetric cells on glass
- Reference measurements
 - ELISA as golden standard
- Measurement setup
- Impedance spectroscopy of SPAG coating
- Outlook

Introduction: Healthy aging

- The project which supports and consolidates medical research, technological development and innovation of chronic geriatric disorders.
- The goal:
 - To decrease costs in the healthcare sector
 - Development of a technological platform for sharing knowledge in the field of healthy aging
- Workpackage T1: Biomarkers
 - Development of a high throughput screening device



Introduction: Sensor principle



- | | | | |
|---|--|---|----------|
|  | Polystyrene |  | Blocking |
|  | Interdigitated electrodes |  | Antigen |
|  | PET (Polyethylene terephthalate)
Or glass |  | Antibody |

Characterization: Contact angle measurements

- Measuring the angle between the liquid–solid interface and the liquid–vapor interface.
- The angle decreases when a substrate is more hydrophilic.
- The angle is a number for the amount of antigen bound.




Characterization:


Antigen coating on substrates

- Polydimethylsiloxane(PDMS) substrate
 - Silicone typically used for microfluidics
 - Hydrophobic
 - Biocompatible
 - Difficult to apply in thin layers
- Polystyrene substrate
 - Same polymer as microplates(MICROLON® 200 96W)
 - Easy to spincoat in layers of $\pm 200\text{nm}$
- Glass
 - Typically the bottom substrate of microfluidics
 - Bad measurements with golden standard




Characterization: Antigen coating on substrates




Characterization: Antigen coating on substrates





Sensor development: Lab-on-card principle





Sensor development: Lab-on-card with impedimetric cell




Sensor development: Lab-on-card with microplate format



Sensor development: Multichannel impedimetric cells on glass




Reference measurement: ELISA as golden standard







Measurement setup


Electrical impedance
spectrometer




Measurement setup: Sensor clamp







Electrical impedance spectroscopy of SPAG coating



Electrical impedance spectroscopy of SPAG coating



Outlook




Outlook: Fully automated immunoassay

Temperature controller



Liquid handler



Electrical impedance spectrometer



Conclusion

- Characterization proves stable coating on polystyrene layer
- Golden standard(ELISA) proves the sensitivity of a measurement in a card
- Full ELISA protocol tests in a card are successful
- Impedance measurements give stable spectra
 - Performing characterization on the interdigitated electrodes, e.g. measuring Silica particles

End of the presentation

Thank you for your attention!