

# Plasma metabolites as first-line responders in lung cancer

Plasma metabolite biomarkers for improved lung cancer diagnosis in patients with solitary pulmonary nodules

Jill Meynen<sup>1</sup>; Peter Adriaenssens<sup>2</sup>; Maarten Criel<sup>1,3</sup>; Karolien Vanhove<sup>1,4,5</sup>; Wouter Marchal<sup>6</sup>; Liesbet Mesotten<sup>1,7</sup>; Elien Derveaux<sup>2</sup>

1 Faculty of Medicine and Life Sciences, Hasselt University, Hasselt

2 Applied and Analytical Chemistry, NMR Group, Institute for Materials Research (imo-imomec), Hasselt University, Hasselt

3 Department of Respiratory Medicine, Ziekenhuis Oost-Limburg, Genk


4 Department of Respiratory Medicine, University Hospital Leuven, Leuven

5 Department of Respiratory Medicine, Algemeen Ziekenhuis Vesalius, Tongeren

6 Applied and Analytical Chemistry, Institute for Materials Research (imo-imomec), Hasselt University, Hasselt


7 Department of Nuclear Medicine, Ziekenhuis Oost-Limburg, Genk

## Problem



PET-CT is an important technique in lung cancer detection and staging. Besides the fact that almost every lesion is detected, this technique often **fails to discriminate benign from malignant PET-positive solitary pulmonary nodules (SPNs)**. This results in prolonged follow-up, delayed diagnosis, and even unnecessary surgery in at least 10% of the patients.


## Solution



This research therefore aims to use **plasma metabolite biomarkers to discriminate between benign and malignant PET-positive SPNs**. The results on this poster will focus on the development of a **RP-HPLC method to detect plasma amino acids**.


## Methods

### I. Plasma collection




Fasted blood samples are collected from all participants and plasma is isolated for further analyses. Samples are stored at -80 °C by UBiLim.

### II. <sup>1</sup>H-NMR & HPLC analysis



Plasma samples are analysed with <sup>1</sup>H-NMR and HPLC. Two HPLC methods will be developed: one to detect the amino acids, one to detect the organic acids.


### III. Statistics



Univariate statistics (unpaired T-test or Mann-Whitney U test) will be performed for all metabolites in the HPLC chromatogram. Multivariate statistics (OPLS-DA) will be performed for all metabolites in the <sup>1</sup>H-NMR spectrum.

## Results

### HPLC method amino acid analysis



1: Aspartic acid; 2: Glutamate; 3: 3-OH-proline; 4: Asparagine; 5: Serine; 6: Glutamine; 7: Glycine; 8: Histidine; 9: Arginine; 10: Threonine; 11: Taurine; 12: Alanine; 13: 1-Methylhistidine; 14: 3-Methylhistidine; 15: Proline; 16: AMQ; 17: Tyrosine; 18: Cysteine & Cystine; 19: Valine; 20: Norvaline (IS); 21: Methionine; 22: Leucine; 23: Ornithine; 24: Lysine; 25: Isoleucine; 26: Phenylalanine; 27: Tryptophan. \*Represent non-identified signals.

## Conclusion

The HPLC method developed to analyse amino acids in plasma represents a crucial step towards identifying potential **plasma biomarkers**. While PET-CT remains invaluable in lung cancer detection, its inability to differentiate **between benign and malignant SPNs** creates a need for complementary diagnostic tools. By combining plasma biomarker analysis with PET-CT, we aim to improve diagnostic accuracy, reduce unnecessary surgery, and improve patient outcomes. Future studies will focus on the validation and identification of specific biomarkers to further refine this approach.

UBiLim  
University Biobank Limburg


LCRC

Ziekenhuis  
Oost-Limburg



Kom op  
tegen  
Kanker

UHASSELT



jill.meynen@uhasselt.be

